

## New solutions for indications with limited interdental space



### The strongest small diameter implant<sup>1</sup>

The Straumann® Standard Plus Narrow Neck Implant (endosteal Ø 3.3 mm / prosthetic Ø 3.5 mm) is the proven benchmark for small diameter implants in the market. It can now be used with confidence for a broader range of indications:

- Single-tooth replacement
- Fixed partial dentures
- Hybrid dentures (retentive anchor available as custom-made product)

### Cutting healing time in half

The Narrow Neck implant is now also available with SLActive, the first chemically active implant surface with hydrophilic properties. SLActive significantly accelerates the osseointegration process. The new Straumann surface offers more security in every indication and improved treatment predictability.

### Easy fabrication of esthetic provisionals

With the new Narrow Neck Temporary Coping customers are now able to time-savily fabricate esthetic provisionals for Narrow Neck implants chair-side. The polymer material is easy to modify and provides a versatile base for crowns and cement-retained bridges, which is ideal for soft tissue conditioning.

Please check the availability of the products with your local Straumann partner.

<sup>1</sup>Comparative fatigue strength measurements of small diameter implants following ISO 14801.